Endurity MRITM

Single-Chamber Pacemaker



Product Highlights - Pacemaker

- Allows patients to undergo 1.5 T or 3 T MRI scans when used with MRI Ready leads from Abbott*
- Physician preferred size and physiologic shape minimize pocket size
- Outstanding longevity provides 14,4 years of service life,¹ which is supported by a 10-year warranty²
- AutoCapture[™] pacing system offers the maximum in threshold adaptability and patient safety with ventricular Beat-by-Beat[™] capture confirmation. The AutoCapture pacing system automatically delivers a 5,0 V backup safety pulse when noncapture is detected, and it may be programmed to either a bipolar or unipolar configuration
- A suite of state-of-the-art features—such as automaticity, Ventricular AutoCapture[™] pacing system and SenseAbility[™] technology—are designed to deliver optimal therapy for patients at implant and throughout their lives
- Real-time electrogram (EGM) waveform, as well as the associated event markers that precede and follow a specific triggering event, can be programmed to automatically record up to 14 minutes of stored EGMs when encountering one or more programmable trigger options
- · 6-month ERI-EOL interval
- An optional, easy-to-use hand-held device
 (SJM MRI Activator™ device) can be used to program
 the device to pre-approved MRI settings pre- and post MRI scan, decreasing the number of workflow steps and
 increasing clinic efficiency

*See MRI Conditional Parameters

Ordering Information

Contents: MRI Ready Pacing System

MODEL NUMBER	DESCRIPTION	DIMENSIONS (H X W X T, MM)	WEIGHT (G)	VOLUME (CC)	CONNECTOR
PM1172	Endurity MRI™ Pacemaker	41 x 50 x 6	19	9.7 (± 0,5)	IS-1

Indications: Implantation is indicated in one or more of the following permanent conditions: syncope, presyncope, fatigue, disorientation due to arrhythmia/bradycardia, or any combination of those symptoms. Rate-Modulated Pacing is indicated for patients with chronotropic incompetence, and for those who would benefit from increased stimulation rates concurrent with physical activity. Atrial Pacing is indicated for patients with sinus node dysfunction and normal AV and intraventricular conduction systems. Ventricular Pacing is indicated for patients with significant bradycardia and normal sinus rhythm with only rare episodes of A-V block or sinus arrest, chronic atrial fibrillation, severe physical disability.

Contraindications: Single-chamber pulse generators are contraindicated in patients with an implanted cardioverter-defibrillator. Rate-Adaptive Pacing may be inappropriate for patients who experience angina or other symptoms of myocardial dysfunction at higher sensor-driven rates. An appropriate Maximum Sensor Rate should be selected based on assessment of the highest stimulation rate tolerated by the patient. Single-Chamber Ventricular Demand Pacing is relatively contraindicated in patients who have demonstrated pacemaker syndrome, have retrograde VA conduction, or suffer a drop in arterial blood pressure with the onset of ventricular pacing. Single-Chamber Atrial Pacing is relatively contraindicated in patients who have demonstrated compromise of AV conduction. For specific contraindications associated with individual modes, refer to the programmer's on-screen help.

Potential Adverse Events: The following are potential complications associated with the use of any pacing system: arrhythmia, heart block, thrombosis, threshold elevation, valve damage, pneumothorax, myopotential sensing, vessel damage, air embolism, body rejection phenomena, cardiac tamponade or perforation, formation of fibrotic tissue/local tissue reaction, inability to interrogate or program a device because of programmer malfunction, infection, interruption of desired device function due to electrical interference, loss of desired pacing and/or sensing due to lead displacement, body reaction at electrode interface, or lead malfunction (fracture or damage to insulation), loss of normal device function due to battery failure or component malfunction, device migration, pocket erosion, or hematoma, pectoral muscle stimulation, phrenic nerve or diaphragmatic stimulation. The following, in addition to the above, are potential complications associated with the use of rate-modulated pacing systems: inappropriate, rapid pacing rates due to sensor failure or to the detection of signals other than patient activity, loss of activity-response due to sensor failure, palpitations with high-rate pacing.

Refer to the User's Manual for detailed indications, contraindications, warnings, precautions and potential adverse events.

Endurity MRI™

Single-Chamber Pacemaker

Product Specifications

PHYSICAL SPECIFICATIONS

Model	PM1172	
Telemetry	Inductive	
Dimensions (mm)	41 x 50 x 6	
Weight (g)	19	
Volume (cc)	9.7	
Connector	IS-1	

Compatible with Merlin@home™ Transmitter

PARAMETER	SETTINGS
Rate/Timing	
(Fixed) (ms)	125; 160-400 in steps of 30; 440; 470 ³
Base Rate (min-1)	30-130 in steps of 5; 140-170 in steps of 10
Mode	VOO(R); VVI(R); VVT(R); Pacing Off
	AOO(R); AAI(R); AAT(R)
Hysteresis Rate (min-1)	Off; 304-150 in steps of 5
Search Interval (min-1)	Off; 1; 5; 10; 15; 30
Cycle Count	1-16 in steps of 1
Intervention Rate (min-1)	Off; 80-120 in steps of 10; Intrinsic +0; Intrinsic
	+10; Intrinsic +20; Intrinsic +30; Same as Base Rate
Intervention Duration (min)	1-10 in 1 minute intervals
Recovery Time	Fast; Medium; Slow; Very Slow
Rest Rate (min-1)	Off; 30-150; in steps of 5
Rate Responsive VREF	Off; Low; Medium; High
Shortest VREF	125-475 in steps of 25

Output/Sensing

ACap™ Confirm⁵ On; Off; Monitor Primary Pulse Configuration Bipolar Bipolar 5,0⁴ Backup Pulse Configuration Backup Pulse Amplitude (V) Search Interval (hours) 8; 24

A or V Pulse Amplitude (V) 0,25–4,0 in steps of 0,25; 4,5–7,5 in steps of 0,5 0,05; 0,1–1,5 in steps of 0,1

A or V Pulse Width (ms)

Unipolar (tip-case); Bipolar (tip-ring) Unipolar Tip (tip-case); Bipolar (tip-ring); Unipolar Ring A or V Pulse Configuration

A or V Sense Configuration (ring-case)

Atrial Sensitivity (mV)

(1.15, 2.15) (1. V Sensitivity (mV)

Ventricular AutoCapture™ Pacing

System

Primary Pulse Configuration Backup Pulse Configuration Backup Pulse Amplitude (V) Search Interval (hours)

SenseAbility™ Technology

Unipolar; Bipolar Unipolar; Bipolar 5,0⁸ 8; 24

(Automatic Sensitivity Control adjustment for atrial or ventricular events) 0,2–1,0 in steps of 0,1

A Max Sensitivity (mV) 0,2-2,0 in steps of 0,1

V Max Sensitivity (mV) Threshold Start (Atrial and Ventricular Post-Sense) 50: 62.5: 75: 100%

(Atrial Post-Pace) 0,2-3,0 in steps of 0,1 mV (Ventricular Post-Pace) Auto; 0,2-3,0 in steps of 0,1 mV (Atrial and Ventricular Post-Sense) 0; 30; 60; 95; 125; 160;

(Atrial Post-Pace) 0; 30; 60; 95; 125; 160; 190; 220 (Ventricular Post-Pace) Auto; 0; 30; 60; 95; 125; 160; 190; 220

AF Management

Decay Delay (ms)

AF Suppression™ Algorithm Lower Rate Overdrive (min-1) 10^{4} Upper Rate Overdrive (min-1) No. of Overdrive Pacing Cycles

Rate Recovery (ms)

Maximum AF Suppression Rate

Atrial Tachycardia Detection Rate (min⁻¹)

Off: On

15-40 in steps of 5

80–150 in steps of 5; 160–180 in steps of 10

110–200 in steps of 10; 225–300 in steps of 25 $\,$

Rate-Modulated Parameters

Maximum Sensor Rate (min-1)	80–150 in steps of 5; 160-180 in steps of 10		
Reaction Time	Very Fast; Fast; Medium; Slow		
Recovery Time	Fast; Medium; Slow; Very Slow		
Sensor	On; Off; Passive		
Slope	Auto (-1); Auto (+0); Auto (+1); Auto (+2); Auto (+3); 1-16 in steps of 1		
Threshold	Auto (-0,5); Auto (+0,0); Auto (+0,5); Auto (+1,0); Auto (+1,5); Auto (+2,0); 1-7 in steps of 0,5		

Stored Electrograms

Off; Low; High		
1; 2; 3		
Off; Low; High		
Off; Low; High		
125-300 in steps of 25		
2; 3; 4; 5; 10; 15; 20		
Off; Low; High		
Off; Low; High		

ar Rate can alternately be High Atrial Rate; they use the same sub-para

Other

V Lead Monitoring	Monitor; Auto Polarity Switch		
$VLow$ Impedance Limit (Ω)	100-500 in steps of 25		
V High Impedance Limit (Ω)	750-2500 in steps of 250; 3000		
Atrial limits apply when implanted in the	e atrium.		
Lead Type Uncoded; Unipolar; Bipolar			
Magnet Response	Off; Battery Test		
NIPS Options			
Stimulation Chamber	Atrial; Ventricular		
Coupling Interval (ms)	100-800 in steps of 101		
S1 Count	2-25 in steps of 1		
S19; S2; S3 and S4 Cycle (ms)	Off; 100-800 in steps of 10 (Fixed or Adaptive)		
Diagnostic Trends	AT/AF Activity, Exercise; Lead Impedance; R (or P) Wave; V (or A) Threshold		

MRI Settings

g-	
MRI Mode	AOO; VOO; DOO; Pacing Off
MRI Base Rate	30-120 min ⁻¹
MRI Pulse Configuration	Bipolar
MRI Pulse Amplitude	5,0 V; 7,5 V
MRI Pulse Width	1,0 ms
MRI Paced AV Delay	25-120 ms

MRI Scan Parameters**

MRI Ready Lead	Lead Lengths	Magnet (Tesla)	Scanner Mode	Scan Region
Tendril™ 2088TC Lead	46, 52, 58 cm		Normal	Full Body
IsoFlex™ Optim™ 1944 Lead	46, 52 cm	1.5 T, 3 T Operating	Operating	
IsoFlex™ Optim™ 1948 Lead	52, 58 cm		Mode	

^{**}Refer to the MRI Ready Systems Manual for more detailed information

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

26020-SJM-END-0914-0004(7) Item approved for international use only.



^{1.} A,V = 2,5 V @ 0,4 ms; 500 ohms; 100% VVI pacing @ 60 bpm; AutoCapture™ Pacing System OFF;

SECMS ON.
Terms and conditions apply; refer to the warranty for details.
Programming options dependent on pacing mode.
The highest available setting for hysteresis rate will be 5 min-1 below the programmed base rate.
Atrial Implants Only.
Values 0,1–0,4 not available in a unipolar sense configuration.
Sensitivity in with senserting a 00 the programmed to the programmed

Values 0,1-0,3 not a variable in a unipola sense configuration. Sensitivity is with respect to a 20 ms haversine test signal. This parameter is not programmable. S1 Burst Cycle is applied at the preprogrammed S1 cycle length.

Customer Support: 46-8-474-4756

[™] Indicates a trademark of the Abbott group of companies.

^{© 2018} Abbott. All Rights Reserved.