

EVOLUTIONARY **WASHER** REVOLUTIONARY **PERFORMANCE**



AMSCO® 7052HP
and **7053HP**
Washer/Disinfectors



Introducing a
New Era in Processing



THE ONLY WASHER/DISINFECTORS WITH A VALIDATED CYCLE FOR da Vinci® INSTRUMENTS*.

AMSCO 7052HP and 7053HP Washer/Disinfectors deliver a validated and cleared process for cleaning da Vinci instruments. The new **RAS** (Robotic Assisted Surgery) Cycle, **RAS 12 Racks** and **Prolystica® Ultra Concentrate HP Enzymatic and Neutral Chemistries** are specifically designed to automate the cleaning and disinfection process — delivering optimal outcomes and saving you time.



Processing for these EndoWrist® Instruments**

- da Vinci Xi®
- da Vinci X®
- da Vinci Si®
- da Vinci S®




TRANSFORMING YOUR PROCESS. FEWER STEPS. MORE PRODUCTIVITY.



with **RAS**

- **Drive consistent outcomes** when using a validated process
- **Reduce time** and variability with automation
- **Improve department productivity** with up to 25% reduction in manual cleaning time per set of 4 instruments

OVERALL DIRECT LABOR SAVINGS
66 minutes*
 using the RAS Cycle to reduce manual cleaning time



DOCUMENTED + RELIABLE + REPEATABLE

= COMPLIANCE


* The RAS Racks are used in the RAS Cycle of the AMSCO 7052HP Single-Chamber Washer/Disinfectors and the AMSCO 7053HP Single-Chamber Washer/Disinfectors for the effective cleaning, rinsing, intermediate level disinfection and drying of reusable da Vinci® X/Xi and S/Si EndoWrist® instruments.

** da Vinci®, EndoWrist® are registered trademarks of Intuitive Surgical

* per set of 4 instruments

PRODUCTIVITY AND THROUGHPUT SECOND TO NONE

AMSCO® 7052HP and 7053HP Washer/Disinfectors take productivity to elevated heights, offering **the industry's highest throughput****** to meet the demands of the busiest SPDs.




DIN Trays per 8 hour day



Continuous throughput

using Cool Touch handles to remove loads immediately after cycle completion



Quicker load preparation

with the versatility of 360° loading

Flexibility to process a wide range of surgical instruments with **multiple rack configurations**



Unlock the **industry's fastest instrument cycle time** with the power of Prolystica® Ultra Concentrate HP Chemistries



Deliver Maximum Cleaning Efficacy with precise chemistry dosing from integrated flow meters


****When compared to current healthcare washer/disinfector models

TECHNICALLY ADVANCED TECHNICIAN FRIENDLY

Stay on Schedule & Informed with cycle progress notification lighting




Easy Cycle Selection
with icon based control panel



Optimal Cleaning
with Complete
Rack Coverage

- Ensure full coverage with cross pattern spray arm with 120° range
- High and low impingement options



The Simplicity of Connectivity


SPM supports best practices by providing alerts to track racks as they move through the sterile processing department. This comprehensive workflow management software receives information directly from the washer and informs staff of required maintenance through on-screen notifications.



SAVING RESOURCES SUSTAINABLE SOLUTIONS

Combining fast, effective cycle times with reduced water consumption minimizes environmental impact while reducing utilities usage.

AMSCO® 7052HP Washer/
Disinfecter uses up to
50% less
water per instrument cycle**



AMSCO® 7053HP Washer/
Disinfecter uses up to
40% less
water per instrument cycle**



**Sources: Performance Specifications relative to STERIS and competitive brands

up to **6** gallons **3** minutes **TOTAL SAVINGS PER CYCLE!**



THE FIRST AND ONLY OF ITS KIND



7052HP

Overall exterior dimensions:	42 x 80-3/4 x 32" (1067 x 2051 x 813 mm)
Load Height:	31" (787 mm) from the floor

7053HP

Overall exterior dimensions:	42 x 80-3/4 x 38" (1067 x 2051 x 965 mm)
Load Height:	31" (787 mm) from the floor

To learn more about the AMSCO 7052HP and 7053HP Washer/Disinfectors, contact your STERIS representative or visit www.STERIS.com/7000HP

 **STERIS**

STERIS
5960 Heisley Road
Mentor, OH 44060-1834 • USA
440-354-2600 • 800-548-4873
www.steris.com