

Library Prep Kit	Index Kit Name	Index Kit Catalog #
TruSeq DNA Nano TruSeq PCR Free	IDT for Illumina - TruSeq DNA UD Indexes v2 (96 Indexes, 96 Samples)	20040870
TruSeq Stranded Total RNA	TruSeq RNA Single Indexes Set A & B (12 Indexes, 48 Samples per set)	A: 20020492 B: 20020493
TruSeq Stranded mRNA	TruSeq RNA CD Index Plate (96 Indexes, 96 Samples)	20019792
	IDT for Illumina - TruSeq RNA UD Indexes v2 (96 Indexes, 96 Samples)	20040871
Illumina Stranded Total RNA Prep, Ligation with Ribo-Zero Plus	Illumina RNA UD Indexes Set A-D, Ligation (96 Indexes, 96 samples per set)	A: 20091655 B: 20091657 C: 20091659 D: 20091661
Illumina Stranded mRNA Prep, Ligation	IDT for Illumina RNA UD Indexes Set A-D, Ligation (96 Indexes, 96 Samples per set)	A: 20040553 B: 20040554 C: 20040555 D: 20040556
Illumina RNA Prep with Enrichment, (L) Tagmentation	Illumina DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20091654 B: 20091656 C: 20091658 D: 20091660
	IDT for Illumina - DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20027213 B: 20027214 C: 20042666 D: 20042667
TruSeq Targeted RNA Expression	TruSeq Targeted RNA Index Kit (48 Indexes, 48 Samples)	RT-401-1001
Illumina DNA PCR-Free Prep, Tagmentation	Illumina DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20091654 B: 20091656 C: 20091658 D: 20091660
	IDT for Illumina - DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20027213 B: 20027214 C: 20042666 D: 20042667

Illumina DNA Prep, (M) Tagmentation <i>[Formerly Nextera DNA Flex]</i>	Illumina DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20091654
		B: 20091656
		C: 20091658
		D: 20091660
	IDT for Illumina - DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20027213
		B: 20027214
C: 20042666		
D: 20042667		
Nextera DNA CD Indexes (96 Indexes, 96 Samples)	20018708	
Illumina DNA Prep with Enrichment, (S) Tagmentation <i>[Formerly Nextera DNA Flex Pre-Enrichment Library Prep and Enrichment Reagents]</i>	Illumina DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20091654
		B: 20091656
		C: 20091658
		D: 20091660
	IDT for Illumina - DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20027213
		B: 20027214
C: 20042666		
D: 20042667		
Nextera XT DNA	Illumina DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20091654
		B: 20091656
		C: 20091658
		D: 20091660
	IDT for Illumina - DNA/RNA UD Indexes Set A-D, Tagmentation (96 Indexes, 96 Samples per set)	A: 20027213
		B: 20027214
		C: 20042666
		D: 20042667
	Nextera XT Index Kit v2 Set A-D (96 Indexes, 384 Samples per set)	A: FC-131-2001
		B: FC-131-2002
C: FC-131-2003		
D: FC-131-2004		
Nextera XT Index Kit (24 Indexes, 96 Samples)	FC-131-1001	
AmpliSeq for Illumina library prep PLUS kits	AmpliSeq UD Indexes for Illumina (24 Indexes, 24 Samples)	20019104
	AmpliSeq CD Indexes Set A-D for Illumina (96 Indexes, 96 Samples per set)	A: 20019105
		B: 20019106
		C: 20019107
		D: 20019167
AmpliSeq CD Indexes Sets A-D for Illumina (384 Indexes, 384 Samples)	20031676	