

# Certified Pre-Owned Equipment

## AMSCO 400 Small Steam Sterilizer



# STERIS CERTIFIED PRE-OWNED AMSCO 400 SMALL STEAM STERILIZER

Our Certified Pre-Owned (CPO) STERIS AMSCO 400 Sterilizer refurbished solution undergoes a precise refurbishing process. STERIS Certified Pre-Owned technicians thoroughly inspect and refurbish each unit to provide and ensure each AMSCO 400 Sterilizer is certified aesthetically and operationally to original equipment manufacturer (OEM) specifications.

---

## STERIS Certified Pre-Owned Sterilizer Certification Process

### Conditional Replacement/Refurbishment

STERIS has developed a comprehensive process that replaces all components with OEM parts as new with the exception of the framing for the stainless, door and the chamber, which is stripped and conditionally painted if needed. Our process includes but is not limited to the following:

- Re-Insulate Chamber
- Re-Finish all Stainless Steel Panels or Replace if Needed
- Replace Steam Manifold
- Replace Vacuum Pump
- Replace Heat Exchanger Assembly
- Replace Water Inlet
- Replace Jacket Trap
- Replace Safety Valve
- Replace Pressure Regulator
- Replace Exhaust Manifold



- Calibrate Check Valve
- Replace Feed Water Valve
- Temperature & Pressure Calibration
- Check Chamber Level Sensor
- Verify Door Operation
- Check Appearance for Leaks and Damage
- Check Relays and Contactors for Arc Damage
- Replace Chamber Strainer
- Replace Conductivity/ Float Probes
- Replace Steam Baffle

# PERFORMANCE AND VALUE.

## EFFICIENT, COMPLIANT STEAM STERILIZATION BUILT FOR PRODUCTIVE PROCESSING.

STERIS steam sterilizers process even the largest loads of instrument sets with effortless operation and efficiency, allowing you to better serve the OR's increased case loads. Our portfolio of models provides you a choice that best fits your needs.

### IDEAL FOR QUICK TURNAROUND IN THE OR OR CSSD

- Meets current AAMI/AORN guidelines
- Process up to 3 - 25lb trays • Validated 4 + 1 prevacuum immediate use cycle
- Customizable cycle names
- Fully programmable cycles to help meet instrument manufacturer processing recommendations



### ORDER A STERILIZER THAT MEETS YOUR FACILITY'S SPECIFIC NEEDS

#### Healthy for the planet and your bottom line.

- Available in two chamber sizes - 16" or 20"
- Install recessed in the wall or free-standing with cabinet
- Use facility steam or choose a fully integrated steam generator
- 20" available with loading car and transfer carriage

#### Sterilizer Dimensions and Chamber Sizes

AMSCO 400 STEAM STERILIZER	INTERIOR CHAMBER DIMENSIONS	MAX. OUTSIDE DIMENSIONS
16X16 Models	16" x 16" x 26" (406 x 406 x 660mm) 108.79 liters	36 7/8"L x 26"W x 74 1/2"H (937* x 660 x 1891mm)
20x20 Models	20" x 20" x 38" (508 x 508 x 965 mm) 249.03 liters	45 5/8"L x 30"W x 74 1/2"H (1159* x 762 x 1891mm)

# STERIS CERTIFIED PRE-OWNED EQUIPMENT

BACKED BY AN EXTENSIVE NETWORK OF  
SKILLED SERVICE TECHNICIANS.



For more information on the STERIS Certified Pre-Owned  
AMSCO 400 Small Steam Sterilizer contact STERIS today.

[healthcare.steris.com](https://healthcare.steris.com)

M13757EN.2025-01, Rev.A  
© 2025 STERIS. All Rights Reserved.

 STERIS