

HOLOGIC®

Introducing the Next Generation of **NovaSure®**

NovaSure V5

Endometrial Ablation



TRUSTED



INNOVATIVE



**PATIENT-CENTRIC
DESIGN**



TRUSTED

Over 20 years of clinical data
97% patient satisfaction²
75% amenorrhea rate at 5 years³
86% of patients avoided hysterectomy at 10 years⁴




INNOVATIVE

Over 3 million procedures performed worldwide has driven continuous design improvements over the years based on user feedback.



PATIENT-CENTRIC DESIGN

To treat a wide range of unique cervical and uterine anatomy through enhanced features.



SureClear™ Technology

Hologic's unique SureClear fluid removal system provides integrated suction through the array by constant tissue contact while simultaneously removing vapour, blood and desiccated endometrium.

Endoform™ Sealing Technology

NovaSure® V5 cervical seal integrates Endoform technology which is designed to accommodate a range of cervical canals and anatomical variability.

AccuSheath™ Markings

NovaSure V5 features a 6mm outer diameter sheath with AccuSheath markings designed to increase accuracy and confidence in device seating and placement.

Bi-Polar RF Array

Quick and controlled bipolar radiofrequency energy is customised to the patient's uterine cavity measurements.



SmartDepth™ technology monitors and measures tissue impedance for a customised depth of ablation for each patient.

ORDERING INFORMATION NOVASURE V5

PRODUCT NUMBER PRODUCT DESCRIPTION

NSV5-003	NovaSure V5 Endometrial Ablation Disposable Device (3-Pack)
NSV5-001	NovaSure V5 Endometrial Ablation Disposable Device (1-Pack)

NovaSure® V5

Endometrial Ablation



CE 2797 EC REP Hologic BV, Da Vincilaan 5, 1930 Zaventem, Belgium

REFERENCES

- ¹ Based on commercial data from Hologic Inc., units shipped from 2004-2018
- ² Baskett TF, Clough H, Scott TA. NovaSure bipolar radiofrequency endometrial ablation: report of 200 cases. Journal of Obstetricians and Gynecologists Canada. 2005;27(5):473-476.
- ³ Gallinat A. An impedance-controlled system for endometrial ablation: Five-year follow up on 107 patients. J Reprod Med. 2007; 52:467-472 1
- ⁴ Herman M, Penninx J, Mol B, Bongers M. Ten-year follow-up of a randomised controlled trial comparing bipolar endometrial ablation with balloon ablation for heavy menstrual bleeding. BJOG 2013;120:966-970

IMPORTANT SAFETY INFORMATION NovaSure® endometrial ablation is for premenopausal women with heavy periods due to benign causes who are finished childbearing. Pregnancy following the NovaSure procedure can be dangerous. The NovaSure procedure is not for those who have or suspect uterine cancer; have an active genital, urinary or pelvic infection; or an IUD. NovaSure endometrial ablation is not a sterilisation procedure. Rare but serious risks include, but are not limited to, thermal injury, perforation and infection. Temporary side effects may include cramping, nausea, vomiting, discharge and spotting. Inform patients to contact you if they experience a possible side effect related to use of this product. If you or someone you know has possibly experienced a side effect when using our product, please contact your physician.

SS-01364-EUR-EN Rev. 002 ©2022 Hologic, Inc. All rights reserved. Specifications are subject to change without prior notice. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative. Hologic, AccuSheath, EndoForm, NovaSure, SmartDepth, SureClear and associated logos are trademarks or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and trade shows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative. For more information contact your local Hologic representative or write to euinfo@hologic.com